




WILFORD

TWO TONE VELVET UPHOLSTERY

ESSENTIA
SERIES



WILFORD | 103



WILFORD | 104



WILFORD | 105



WILFORD | 106




WILFORD | 107



WILFORD | 108



WILFORD | 109



WILFORD | 110



WILFORD | 111



WILFORD | 112

"Note that actual product colors may vary slightly in digital view due to lighting and screen settings."



WILFORD | 115



WILFORD | 114



WILFORD | 116



WILFORD | 117



WILFORD | 113



WILFORD | 118



WILFORD | 121



WILFORD | 120



WILFORD | 122



WILFORD | 123



WILFORD | 124



WILFORD | 125

"Note that actual product colors may vary slightly in digital view due to lighting and screen settings."



WILFORD | 128



WILFORD | 126



WILFORD | 129



WILFORD | 127



WILFORD | 134



WILFORD | 133



WILFORD | 130



WILFORD | 131



WILFORD | 132



WILFORD | 137



WILFORD | 136



WILFORD | 135

"Note that actual product colors may vary slightly in digital view due to lighting and screen settings."



WILFORD | 113

WILFORD

TWO TONE VELVET UPHOLSTERY



COLLECTION	WIDTH	GSM	COMPOSITION	MARTINDALE	COLOUR FASTNESS TO LIGHT	WASH CARE	END USE
WILFORD	140	250	100% Polyester	40000	4		
WASH CARE :- Wash at 40 degree Tumble dry low heat Only non-chlorine bleach Dry clean Don't Iron				END USE :- Upholstery Complimentary Sample Not For Sale Color May Vary From Batch To Batch			

Wash care instructions need to be followed for prolonged usage of fabric, otherwise will damage the fabric. Test results were gathered from various samples and therefore represent average indicative values. Discolouration caused by clothing fabrics such as jeans does not representa material defect.Minor colour deviations & discolouring in context of permissible industrial tolerance standards.



WILFORD

TWO TONE VELVET UPHOLSTERY

ESSENTIA
SERIES